

CLASSIC

Therapeutic graduated compression stockings crafted with superior-quality yarns offer a line of reliable, durable garments with fine knitwork and meticulous attention to detail.

Standard length knee-high open-toe

K2 compression stockings

23-32 mmHg

Composition:

70% PA - 30% EA

PRODUCT CODE:

M2450A



Available colours:

● Beige

Short length knee-high open-toe

K2 compression stockings

23-32 mmHg

Composition:

70% PA - 30% EA

PRODUCT CODE:

M2450AS



Available colours:

● Beige

Thigh-high open-toe K2 compression

stockings with hypoallergenic

silicone top bands

23-32 mmHg

Composition:

70% PA - 30% EA

PRODUCT CODE:

M2470A



Available colours:

● Beige

Open-toe K2 compression

slimming

tights

23-32 mmHg

Composition:

70% PA - 30% EA

PRODUCT CODE:

M2480A



Available colours:

● Beige

SIZE TABLE		1 - S	2 - M	3 - L	4 - XL	5 - XXL
Ankle circumf.	cm	19 - 21	21 - 23	23 - 25	25 - 27	27 - 29
Calf circumf.	cm	29 - 34	32 - 37	35 - 40	38 - 43	42 - 47
Thigh circumf.	cm	39 - 44	43 - 48	47 - 52	51 - 56	55 - 60
LENGTH STANDARD SIZES						
From ground to knee	cm	38 - 39	39 - 40	39 - 40	40 - 41	40 - 41
From ground to gluteus	cm	71 - 72	72 - 73	73 - 74	74 - 75	75 - 76
LENGTH SHORT SIZES						
From ground to knee	cm	35 - 36	35 - 36	36 - 37	36 - 37	36 - 37
From ground to gluteus	cm	65 - 66	66 - 67	66 - 67	67 - 68	67 - 68
FIBER CODE CAPTION		WASH			CERTIFICATES	
CO	Cotton					
EA	Elastane					
ME	Metallic (Silver)					
PA	Nylon					