

ELBI PREMIUM

Epicondylitis brace SRT 207

Indications

Used for epicondylitis, tendonitis, elbow joint overstressing (tennis elbow / golfer's elbow).

Description

Brace made of textile, comprising of:

- Targeted compression pad at the epicondyle to absorb forces exerted on the tendon
- Adjustable strap with buckle. Hook and loop fastener.

Width: 8 cm. Due to the adjustable strap this orthosis can be tightened to achieve the optimal compression over the tendon.

Composition

Neoprene-coated polyester fabric.



Sizing: One size (adjustable)





 compression pad at the epicondyle, that absorb forces on the tendon.



 hook and loop fastener with buckle