




PRIMO 1 Leather Orthotic Insole

These light, comfortable and breathable leather orthotic insoles are designed to treat and offer relief from a wide range of common foot pain problems. Reduces metatarsal pressure, and reduce strain on Achilles' tendon.

EXCELLENT ARCH SUPPORT:



ORDER CODE:PRIMO 1

- Tanned leather top cover, breathable and hygienic
- Latex rubber heel pad offers excellent shock absorption and cushioning
- Metatarsal pad-supports the metatarsal arch and helps redistribute weight from the ball of the foot
- High arch support helps to reduce overpronation
- Supports a high arch
- Redistributes foot pressures
- Reduces strain on Achilles' tendon
- Improves posture and stability

SIZE RANGE:

EU SIZE RANGE	35-36	37-38	39-40	41-42	43-44	45-46
UK SIZE RANGE	2.5-3.5	4-5	6-6.5	7-8	9-10	10.5-11

KEY FEATURES:

- Wipe Clean with a Damp Cloth**
- Arch Support**
- Heel Support**
- Shock Absorption**
- Sweat-Absorption & Deodorisation**