

PROGENIE Children's Orthotic Insole

These comfortable orthotic insoles for children are designed to treat and offer relief from a wide range of foot pain problems common in children.


EXCELLENT SUPPORT:

3 qtr. Length

3 qtr. Length to allow fitting into various shoe sizes

Side Support

Moulded side support to prevent rolling



Arch Support

Excellent high arch support

Stable Heel Post

Neutral heel post gives a stable base of support

Heel Support

Deep heel cup to stabilise the heel

SIZE RANGE:

SIZE	XXS	XS	S	M	L	XL
EU SIZE RANGE	C 6 - C7	C 7.5 - C8.5	C9 - C10	C10.5 - C11	C12 - C13	C13.5 - UK1.5

KEY FEATURES:



Wipe Clean with a Damp Cloth



Arch Support



Heel Support

ORDER CODE: PROGENIE

- Moulded side support increases stability by cradling the foot and preventing foot roll
- 3 qtr. length is easier to swap between shoes
- Excellent arch support to help improve posture
- neutral heel post - stable base of support
- Deep heel cup to stabilise the heel