


PURA 5 Orthotic Diabetic Insole

These light, comfortable and breathable orthotic diabetic insoles are designed to treat and offer relief from a wide range of common foot pain problems including mild diabetic neuropathy and rheumatoid arthritis of the foot and ankle.

LAYER STRUCTURE:



ORDER CODE: PURA5

- Gentle TPU arch reinforcement reduces overpronation & supports a high arch
- Redistributes foot pressures
- Improves posture and stability
- Reduces metatarsal pressure
- Reduces strain on Achilles' tendon
- Gel anti-shock heel pad - cushions and reduces shock of impact on foot and ankle
- The self moulding Plastizote® top cover for sensitive feet also reduces friction
- Provides soft cushioning
- Discourages bacterial growth

*Please note that if you have neuropathy, you should consult with your doctor or other relevant healthcare professional before starting to wear any insole. You may not be able to feel whether your insoles are uncomfortable or if they are causing damage. It's important therefore that you check your feet daily for any signs of damages or change.

SIZE RANGE:

SIZE	S	M	L
EU SIZE RANGE	35 - 39	37 - 41	42 - 47

KEY FEATURES:



Wipe Clean with a Damp Cloth



Arch Support



Heel Support



Shock Absorption



Sweat-Absorption & Deodorisation