

MANU AKTIV PLUS

Mobile wrist brace

SRT 202

Indications

Used for sprains, tenosynovitis, intense sporting activities, and minor injuries.

Description


The brace is made of textile, with a buckle fastening strap. Hook and loop fastener. Width: 8 cm.

This orthosis is a support for the radiocarpal joint and can be used both as orthopaedic treatment in minor ailments and for recovery in the post immobilization period.

Due to the fastening band and the buckle, it can be fixed at the optimal fit and can be used during intense sporting activities which overload the radiocarpal joint, such as golf.

Composition

The three-layer laminated textile behaves like a moisture regulator, ensuring a quick transfer of the moisture in an outward direction.



Air Transfer



Water Transfer



- Available in Black



- support for the radiocarpal joint



- strap with buckle fastening

Sizes	1	2	3
Wrist circumference (cm)	15-18	18-21	21-25
Reference	SRT202M1	SRT202M2	SRT202M3



* Measurement taken from the circumference of the wrist.