

MANU DYNAMIC

Mobile wrist brace

SRT 206

Indications

Used for tenosynovitis, sprain of the 1st metacarpal bone, and intense sporting activities.

Description

Brace made of textile, with Velcro fastening. Width: 9,5 cm. The orthosis is a support for the radiocarpal joint and can be used both as orthopaedic treatment in minor injuries and for recovery in the post immobilisation period.

Composition

Neoprene coated polyester fabric.



- support for the wrist joint



- Adjustable strap with buckle fastening for optimal fit



- Available in Black

Sizing:

One size
(adjustable)