

OMBIFIX

Abdominal support for umbilical hernia SRT 110

Indications

Used for treatment of an umbilical hernia, abdominal weakness, or trauma, and for post-op recovery.

Description

The support is made of elastic fabric and provided with one rigid, targeted pressure pad (polyethylene foam, Ø 12cm), on the front at the level of the navel. Lateral closure with hook and loop straps. Width: 16 cm.

It is a traditional binder that due to the material, allows normal thoraco abdominal breathing, while ensuring uniform compression.

The semirigid pad provides protection and compression to the umbilical region. The support can be used for umbilical hernia, both before and after surgery.

Composition

Elastic fabric: 50% polyester, 33% polyamide, 17% elastane, with a ventilated structure that allows high air permeability in both directions. Abdominal pad: Polyethylene foam lined with 100% soft cotton.



Air Transfer




Available in Black or Beige





Semirigid targeted compression pad at the level of the navel



 elastic fabric that conforms to the shape of the body





^{*}Measurement taken from the circumference of the waist.