

STOMATEX DUO

Abdominal stoma support SRT 135

Indications

Used after abdominal surgery, for people with a colostomy or urostomy. Used to prevent and treat parastomal hernia.

Description

The support is made of elastic fabric. The frontal part is fitted with an adjustable opening for the colostomy bag.

Hook and loop fastener on the front side.
Width: 20 cm. It is a traditional elastic corset that due to the elastic fabric composition allows normal thoraco abdominal breathing.

The abdominal section has adjustable Velcro closures to customise the size of the opening for the colostomy bag. This optimises fit and therefore can increase support around the stoma and the surrounding abdominal muscles.

Composition

Elastic fabric: 50% polyester, 33% polyamide, 17% elastane, with a structure that allows high air permeability in both directions.




Air Transfer



- Available in Black or Beige




- elastic corset design provides comfortable and snug compression whilst not impeding function
- the bag can be easily attached and accessed



- soft elastic fabric for maximum comfort
- adjustable opening for the bag

Sizes	1	2	3	4	5	6
Waist circumference (cm)	65-75	75-85	85-95	95-105	105-115	115-125
Reference	SRT135M1	SRT135M2	SRT135M3	SRT135M4	SRT135M5	SRT135M6



*Measurement taken from the circumference of the waist.