

STOMATEX RING

Abdominal stoma support

SRT 136

Indications

Used after abdominal surgery for people with colostomy or urostomy. Used to prevent and treat parastomal hernia.

Description

The support is made of elastic fabric. The frontal, rigid part is fitted with an opening for the colostomy bag, made with a 25 mm wide and 75 mm inside diameter polyethylene ring (PE1).

Hook and loop fastener on the front side. Width: 20 cm. It is a traditional elastic corset that due to the elastic fabric composition allows normal thoraco abdominal breathing.

Composition


Elastic fabric: 50% polyester, 33% polyamide, 17% elastane, with a structure that allows high air permeability in both directions.




Air Transfer



- Available in Beige




- elastic corset design provides comfortable and snug compression whilst not impeding function
- the bag can be easily attached and accessed



- soft elastic fabric for maximum comfort

Sizes	1	2	3	4	5	6
Waist circumference (cm)	65-75	75-85	85-95	95-105	105-115	115-125
Reference	SRT136M1	SRT136M2	SRT136M3	SRT136M4	SRT136M5	SRT136M6



*Measurement taken from the circumference of the waist.