

# TRIAFIX HERNIA 02

## Hernia Support

SRT 145

### Indications

Truss for small - medium reducible inguinal hernia. For pre or post-op support.

### Description

Textile support with adjustable waist and leg straps for comfortable and secure fit.

Includes 1 removable compression pad to increase effectiveness of the support.

Hook and loop fastener on the side. The orthosis is discrete and can be comfortably worn under clothing.

### Composition

Breathable polyester textile. Elastic strap: 50% polyester, 33% polyamide, 17% elastane.

Compression pad: Polyethylene foam.



- Adjustable waist and leg strap



- Compression pad inserts into inner pocket



- Available in Black

Sizes	1	2	3
Hip circumference (cm)	75-95	95-115	115-135
Reference	SRT145M1	SRT145M2	SRT145M3



\*Measurement taken from the circumference of the hip

\*Available in Left or Right, please specify when ordering.