

TRIAGEN

Mobile knee brace

SRT 319

Indications

Used for medial and lateral instability of the knee, sprains and strains of the knee, inflammation, arthritis and for postoperative recovery.

Description


Hinged wrap-around brace made of textile, consisting of:

- 2 removable polycentric hinges, one on each side
- 2 straps with hook and loop fastening


Hook and loop fastening on the front side.
Length: 32 cm. The wraparound design of the brace allows adjustability for oedema, swelling and deformity.

Composition

Three-layer laminated moisture wicking textile.



Air Transfer



Water Transfer



- the knee brace fastens at the front for ease of fitting and adjustability
- polycentric metal hinges for moderate medio-lateral stability




- secure hook and loop straps to improve suspension and reduce slipping



• Available in Black

Sizes	1	2	3	4	5	6	7	8
Knee circumference (cm)	33-35	35-37	37-39	39-41	41-43	43-45	45-47	47-49
Reference	SRT319M1	SRT319M2	SRT319M3	SRT319M4	SRT319M5	SRT319M6	SRT319M7	SRT319M8



*Measurement taken from the circumference of the knee is measured at the middle of the knee cap